Dreame Chrona Pro Max Photon Anti-Aging Skincare Set

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WARNING

To avoid electric shock, fire, or other accidents, carefully read this user manual before use and keep it for future reference.

- Before using the Anti-Aging Skincare Set, you must first perform a sensitivity test (refer to the Photosensitivity Test section). If you experience allergic reactions to the Anti-Aging Skincare Set, do not use this product. This product is intended for home healthcare use only and may be operated by the patient themselves.
- Do not use the Anti-Aging Skincare Set in ways not explicitly specified in these instructions. Any use of the Anti-Aging Skincare Set other than those explicitly specified in these instructions has not been tested, and the associated risks are unknown. This product is not designed for use by individuals with physical, sensory, or mental limitations, or those lacking experience and knowledge, unless they are under adequate supervision or have been instructed on safe use to avoid fire or burn risks.
- If you feel the product overheating or your skin becoming uncomfortable due to excessive heat, immediately turn off and remove the Anti-Aging Skincare Set.
- This product operates wirelessly during use. Connect the charging cable only during charging.
- During charging, the charging cable may pose a risk of entanglement or strangulation. Always keep the product and charging cable out of reach of children at all times.
- Always remove makeup thoroughly before using the Anti-Aging Skincare Set for treatment. Certain cosmetic ingredients may cause photosensitivity or reduce treatment effectiveness. Using the Anti-Aging Skincare Set without thoroughly removing makeup may cause skin reactions or reduce treatment effectiveness. This product is not suitable for use by children. Ensure proper supervision to prevent children from playing with the product, to avoid the risk of fire and burns.
- Never use the Anti-Aging Skincare Set near water or in environments where it may be splashed with water. If you suspect the Anti-Aging Skincare Set has been exposed to water, do not use it. Do not subject the product to strong impacts to avoid the risk of lamp damage.
- Keep the Anti-Aging Skincare Set away from hot surfaces. If necessary, contact the manufacturer for assistance with installation, use, or maintenance of the product, or to report abnormal operation or events.
- Do not use the Anti-Aging Skincare Set if the product, remote control, or USB Type-C charging cable is damaged. Do not leave the product unattended when powered on to prevent fire, burns, or other hazards.
- When removing the product from the lowest/highest storage temperature environment, it is recommended to let it sit for 2 hours before use.
- Do not stare directly at the light or allow it to shine into the eyes during use.
- Do not use the Anti-Aging Skincare Set near devices that use wireless technology.

Precautions

- This product contains silicone, plastic components, and Velcro material. Those known to be allergic to any of the above materials should discontinue use and consult a medical professional.
- This product has a memory function. It automatically remembers and restores the last used treatment mode, intensity level, and treatment time each time it is powered on.
- The Anti-Aging Skincare Set cannot be used in charging mode.
- Clean the product with it turned off before each use.
- Do not use any cleaning products other than those mentioned in the cleaning instructions. Using other cleaning products may damage the Anti-Aging Skincare Set.
- Store the Anti-Aging Skincare Set in a dry environment at room temperature, away from direct sunlight.
- When removing the product from the lowest/highest storage temperature environment, it is recommended to let it sit for 2 hours before use.
- Do not stare directly at the light or allow it to shine into the eyes during use.

Electrical Safety


- Do not damage, modify, bend, twist, or heat the USB Type-C charging cable. A damaged charging cable may cause fire, electric shock, or malfunction.
- Do not place or operate the device near water or other liquids. Otherwise, it may result in fire or electric shock.
- This product uses a magnetic detachable battery pack. Before using the product for the first time, fully charge the battery pack.
- Only use the dedicated battery pack provided with the product.
- Do not attempt to open, modify, replace, or use unapproved battery packs. Improper handling may cause overheating, leakage, or explosion.
- Do not dispose of batteries carelessly. Follow local disposal regulations for disposal.
- Only use a DC 5V $\approx$ 1A charger to charge this product (maximum output power 3.5W).
- Do not use chargers of unknown origin, poor quality, or with an output power exceeding 30W, as this may cause product overheating or damage. For maximum safety, it is recommended to use a charger that complies with IEC 60601-1:2005 + AMD2:2020 standard and has the following specifications: Input: 100-240VAC Output: DC 5V $\approx$ 1.0A.

The product body is marked "Do not use while bathing or showering".

Product Introduction

EN

1 Power button

2 Battery indicator

4 green lights on: Battery 100%
(Fully Charged)
3 green lights on: remaining battery
approximately 75%
2 green lights on: remaining battery
approximately 50%
1 green lights on: remaining battery
approximately 25%
Red light flashing: battery below 10%,
low battery

3 Display screen

View selected mode/intensity/timer

4 Timer key

5 minutes/10 minutes/15 minutes


5 Intensity level key

Level I (low intensity)
Level II (medium intensity)
Level III (high intensity)

6 Power/mode
selection key

Five modes: Firming mode,
Brightening mode, Spot-fading
mode, Smoothing mode, Multi-zone
all-in-one mode

7 Type-C port

1 Neck device

2 Neck device battery indicator

4 green lights on: battery 100% (fully charged)
 3 green lights on: remaining battery approximately 75%
 2 green lights on: remaining battery approximately 50%
 1 green lights on: remaining battery approximately 25%
 Red light flashing: battery below 10%, low battery

3 Remote control display

View selected mode/intensity/timer

4 Timer key

5 minutes/10 minutes/15 minutes

5 Intensity level key

Level I (low intensity)
 Level II (medium intensity)
 Level III (high intensity)

6 Power/mode selection key (Remote control)

Neck device: only one care mode available

7 Type-C port

Note: All product images, accessories, and user interface illustrations in this manual are for reference only. Due to product updates and upgrades, the actual product may differ slightly from the illustrations.

How to Use

Setting Up the Device

1. Remove the device and all included accessories.
2. The mask is equipped with two magnetically removable batteries (each with a capacity of 1000mAh, totaling 2000mAh), and the neck cover contains a 1000mAh battery. Before first use, fully charge both batteries and the remote control. For detailed instructions, please refer to the "Charging Instructions" section.

Photosensitivity Test

Follow the steps below to test whether you have a sensitivity reaction to the treatment light:


1. Press the power button to start the device. Use the remote control to turn the device on or off and select the desired treatment mode.
2. Initial test: Place the beauty mask close to a small area of skin on your forearm and complete one full skin treatment session (three minutes) in the selected mode.
3. Switch to another mode and perform a second treatment session (three minutes) on the same skin area.
4. After completing the test, wait for 24 hours. During this period, no adverse reactions should occur in the test area.

Definition of adverse reactions: Pain or discomfort at the treatment site, as well as abnormal phenomena such as persistent skin redness.

Charging Guide


The device is equipped with two detachable 1000mAh battery packs and a wireless remote control (with a built-in 280mAh battery).

Before first use, follow these instructions:

Battery Pack Charging Instructions


- Before charging, please remove the battery packs from the mask.
- Each battery pack takes approximately 2 hours to fully charge.
- Actual charging time may vary depending on the remaining battery level.

Charging Guide

The neck device is equipped with a built-in rechargeable lithium battery (neck device battery capacity: 1000mAh).

Before first use, follow these instructions:

Battery Pack Charging Instructions

- Use the included Type-C charging cable to connect the charging port of the neck device, and plug the other end into a 5V/1A power adapter.
- During charging, the power indicator on the neck device flashes red (breathing effect); when fully charged, it stays solid green.
- Full charging time is approximately 100 minutes; actual charging time may vary depending on the remaining battery level.

Remote Control Charging Instructions

- This wireless remote control has a built-in 280mAh rechargeable battery.
- Full charging time is approximately 2 hours.
- Actual charging time may vary depending on the remaining battery level.

Lithium Battery Maintenance Recommendations

- Avoid charging after the battery is completely depleted.
- Avoid overcharging the battery; do not exceed 8 hours of charging.
- If the device is not used for a long period, it is recommended to charge the battery at least once every three months.

EN Getting Started

Preparation Before Use

1. Clean your face (if wearing makeup, remove it first and wash your face).
2. This product can be used both in the morning and evening. It can be applied after cleansing or integrated into your daily skincare routine.
3. Ensure the battery is fully charged before use.

Steps for Use

1. Press the power button lightly to turn on the mask (it can be turned on when the green light flashes).
2. After adjusting the headband to a suitable tightness, wear the mask so that the inner side fits snugly against the face.

Headband adjustment guide

1. When using in a seated position, adjust the headband to minimize pressure on the eye sockets and ensure the silicone eye pads fit closely around them to secure the product.
2. If you feel excessive pressure from the silicone eye pads on your eye sockets, loosen the headband straps on both sides appropriately.
3. If there are gaps between the eye pads and eye sockets causing light leakage, tighten the headband straps on both sides.
3. If the mask easily slips off and cannot maintain its correct position, tighten the headband straps as needed.
5. Use the buttons on the remote control to select the desired care mode.
6. If you feel discomfort during wear, try adjusting the headband or remove it to rest until comfortable.
7. The product will automatically shut off after the set mode time has elapsed.

Preparation Before Use

1. Cleanse the neck skin (if wearing makeup or applying skincare products, it is recommended to cleanse first)
2. This product can be used both in the morning and evening, and can also be integrated into your daily skincare routine.
3. Ensure the battery is fully charged before use.

Steps for Use

Without remote control:

1. Wrap the neck device around your neck, adjust the Velcro elastic strap to a comfortable position, ensuring the lamp beads side fits against the neck skin.
2. Press and hold the power button on the neck device for 2 seconds to turn it on.
3. The default timer is 10 minutes, with high light intensity, initiating neck care.
4. After the set time ends, the neck device will automatically power off.

With remote control, the neck device must be paired with the remote control:

Pairing steps: Press and hold the power button on the LED neck device for 5 seconds to enter pairing mode. The indicator light on the neck device will switch from green to red, then quickly flash between green and red. Lightly press the power button on the remote control to turn it on. A beep sound indicates successful pairing.

1. Wrap the neck device around your neck, adjust the Velcro elastic strap to a comfortable position, ensuring the lamp beads side fits against the neck skin.
2. Before daily use, press and hold the power button on the neck device for 2 seconds to turn it on.
3. Press and hold the power button on the remote control to enter care mode.
4. During care, you can adjust the duration (5/10/15 minutes cycle) and intensity (I/II/III) by pressing the timer button on the remote control as needed.
5. Once the set time ends, the neck device will automatically shut off.

Wearing Adjustment Guide

When wearing, ensure the neck device fits snugly against the neck but not too tight to avoid discomfort.

If you feel excessive pressure, loosen the Velcro strap appropriately. If the neck device tends to slip off, tighten the Velcro strap appropriately.

Precautions

- The remote control should be used within 2 meters of the product.
- Avoid intense physical activity while wearing the product.
- During use, you can lie down or lean on a sofa or bed to relax your body.

Remote Control Operation Instructions

1. Connect the included charging cable to the Type-C port of the remote control.
2. Plug the USB end of the charging cable into a 5V 1A power adapter for charging.
3. Charging status display:
 - "88" flashing → Charging in progress
 - "88" Steady On → Charging Complete
 (Unplug the charging cable before use)
4. Low Battery Alert: When the remote control battery is low, "88" will flash continuously until the screen goes off. Please recharge promptly.

Since the LED mask uses a wireless RF remote control (the remote control sends RF signals, and the main unit receives them), you need to pair the remote with the LED mask upon first use.

Pairing Steps: Press and hold the power button on the LED mask for 5 seconds to enter pairing mode. The indicator light on the mask will switch from green to red, then rapidly flash between green and red.

Power On/Off Operation

EN

Mask Mode Selection

Press the Power/Mode Select key (this key serves as both the power switch and the mode selector) to toggle among the following care modes:

Mode 1: Firming Mode Red Light (630nm) + Near-Infrared Light (850nm) + Near-Infrared Light (1070nm)

Effect: Target early signs of aging and skin laxity, with triple light energies working synergistically to penetrate deep into the skin and promote collagen regeneration. Significantly improve elasticity and firmness, accelerate cell repair processes, restoring smoothness, fullness, and youthful radiance to the skin.

The 1070nm wavelength has strong penetration, reaching the deep dermis, subcutaneous tissue, and even the superficial muscle layer. Compared to shallower light waves, it can act on deeper skin structures, helping to improve fine lines, sagging, and enhance skin firmness and elasticity.

Mode 2: Brightening Mode Yellow Light (590nm)

Effect: Soothe and revitalize the skin, especially suitable for recovery after sun exposure or fatigue. Help brighten skin tone, even out texture, and restore a healthy glow.

Mode 3: Spot-fading Mode Green Light (530nm)

Effect: Calm and balance sensitive or stressed skin, reduce redness, and promote an even complexion. Leave the skin in a naturally balanced and refreshed state, especially suitable for reactive or sensitive skin.

Mode 4: Smoothing Mode Blue Light (465nm)

Effect: Commonly used in skincare routines to help soothe visible redness, reduce the appearance of shine, and support a clearer complexion. Maintain a refreshing and balanced skin condition, promoting healthier and more refined skin texture.

Mode 5: Multi-Zone All-in-One Mode (Yellow Light 590nm + Green Light 530nm + Red Light 630nm + Near-Infrared Light 850nm + Near-Infrared Light 1070nm)

Wavelength Combination: Yellow Light (590nm) + Green Light (530nm) + Red Light (630nm) + Near-Infrared Light (850nm + 1070nm)

Effect: Integrate multi-band light energy to synergistically achieve four skincare goals: anti-wrinkle, brightening, spot-fading, and purification. Revitalize the skin base, even out skin tone, reduce pigmentation, and improve overall skin texture. Serve as a daily all-in-one facial care routine with both restorative and preventive benefits, continuously radiating a youthful glow.

Neck Device Treatment Mode Instructions

Mode: Firming Mode

Wavelength Combination: Red Light (630nm) + Near-Infrared Light (850nm) + Near-Infrared Light (1070nm)

Effect: Specifically designed for neck skin care, using the gentle companionship of triple light energy to deeply penetrate and deliver abundant hydrating energy. Help maintain the firmness and elasticity of neck skin, turning daily care into a pleasurable experience. With consistent use, the neck skin becomes increasingly smooth and soft, restoring a delicate and radiant youthful glow.

Special Notes

1. Once the set time ends, the neck device will automatically shut off.
2. Turning off the remote control will also power off the neck device.
3. When the neck device buttons are operated again some time after a treatment mode is selected via the remote control, the device will default to a 10-minute timer. At this point, the timer display on the neck device and the remote control will be inconsistent. Therefore, it is recommended that you select a mode with the remote control and start the light therapy directly.

Cold Compress Guide

You can choose an appropriate cold compress method based on your personal needs and skin condition:

I. Preparing the Cold Compress Pack

Gentle Cold Compress

If you need soothing, relaxation, or fatigue relief, place the cold compress pack flat in a refrigerator (temperature approximately 0°C to 10°C) and take it out after about 15 minutes for use. The temperature of the cold compress pack is approximately 20°C.

*Note: It is recommended to store the cold compress pack in a storage bag inside the refrigerator on a daily basis.

Deep Cold Compress

If you need rapid cooling or enhanced cold compress effects, place the cold compress pack flat in a freezer (temperature approximately -15°C to -20°C) and take it out after about 10 minutes for use. The temperature of the cold compress pack is approximately 15°C.

Storage Tips

It is recommended to store the cold compress pack in the refrigerator on a daily basis, and only place it in the freezer when deep cold compress is needed. Do not place it in the freezer for more than 15 minutes, as this may result in icy slush.

II. Instructions for Use

After removing the cold compress pack from the refrigerator, attach the magnetic side of the pack to the corresponding magnetic side of the mask, then apply it to the desired area.

It is recommended to limit each cold compress session to 15 to 20 minutes, with adjustments based on personal tolerance. If continuous cold compress is required, wait at least 1 to 2 hours before the next session.

III. Function Description of Two Cold Compress Packs

This product set includes two accessories: a full-face cold compress pack and a local cold compress pack, each designed to meet cold compress needs for different areas:

Full-face Cold Compress Pack


Suitable for large-area cold compress on the face, helping to alleviate overall heat sensation, redness, or mild swelling. Ideal for uniform cooling across the entire face.

Tip: The full-face cold compress pack comes with straps. Due to individual differences in facial contours, if the pack does not fully fit the face, it can be secured with the straps.

Local Cold Compress Pack (fit the lower eye area)

Designed specifically for the eye area, fit the lower eye contour, to relieve puffiness, dark circles, or fatigue. Especially suitable for precise cold compress around the eyes after staying up late, prolonged eye strain, or beauty treatments.


Cold Compress Guide

Full-face Cold Compress Pack

Local Cold Compress Pack

After removing the cold compress pack from the refrigerator, attach the magnetic side of the pack to the corresponding magnetic side of the mask, then apply it to the desired area. It is recommended to limit each cold compress session to **15–20 minutes**, with adjustments based on personal tolerance. If continuous cold compress is required, wait at least 1 to 2 hours before the next session.

How to Use

Use with Beauty Mask

After using the LED beauty mask, choose the full-face or local cold compress pack based on facial condition to help soothe and calm the skin.

Standalone Use

The cold compress pack can also be used independently without the mask. It is suitable for daily eye fatigue relief, facial swelling, minor bump cooling, and other scenarios, offering flexibility and convenience.

* If you wish to remove the cold compress pack, you can still use the mask alone for light therapy.

Precautions for Cold Compress Pack

If you experience numbness, tingling, or skin whitening during cold compress, discontinue use immediately. Use with caution or under a doctor's guidance if you have skin damage, eczema, or circulatory disorders.

When children, the elderly, or those with reduced sensitivity use the cold compress pack, adult supervision and close observation are required.

Intensity Level Adjustment

The mask offers the following levels for selection:

Intensity 1 (Low)

Intensity 2 (Medium)

Intensity 3 (High)

Timer Setting

Press the timer setting key to select the following usage durations:

- 5 minutes
- 10 minutes
- 15 minutes

Memory Function Description

When used for the first time, the product is set by default to: Mode 1 (Revitalizing Red Light), Intensity Level 3 (High), Timer 10 minutes. After each use, the product will automatically remember the last used treatment mode, intensity level, and timer setting.

Special Notes

1. When the set treatment time ends, the face mask will automatically turn off.
2. Turning off the remote control will also power off the face mask.
3. When the mask buttons are operated again some time after a treatment mode is selected via the remote control, the mask will default to a 10-minute timer. At this point, the timer display on the mask and the remote control will be inconsistent. Therefore, it is recommended that you select a mode with the remote control and start the light therapy directly.

Cleaning and Maintenance

1. Store the face mask in a normal temperature and dry environment. It is recommended to use the included storage bag.
2. Do not store the product in a dusty, direct sunlight, or humid environment.
3. Do not disassemble or modify internal components on your own.
4. Before each use, turn off the product power and wipe the contact surface of the mask according to the following cleaning instructions.
5. Avoid excessive bending or twisting of the silicone part of the mask.
6. Do not apply cosmetics directly to the surface of the mask.
7. Keep the product away from oils, lotions, and perfumes during use.

Cleaning Instructions

1. The main body of the mask is made of silicone. After each use, gently wipe the surface with a soft, lint-free cloth dampened with clean water.
2. For deep cleaning, use a mild, alcohol-free, and fragrance-free neutral soap diluted with clean water to gently wipe the surface, then wipe clean with a damp cloth.

Daily Cleaning of Cold Compress Pack

1. After each use, gently wipe the surface of the cold compress pack with a soft damp cloth or sponge using a small amount of neutral detergent (e.g., dish soap) to remove sweat or dirt.
2. Thoroughly wipe off any detergent residue with clean water to prevent chemical corrosion of the PVC material.
3. Hang or lay the cold compress pack flat in a cool, ventilated area to air dry naturally. Do not expose to direct sunlight, use a dryer, or blow hot air, as this may accelerate aging.

Recommendations for of Storage and Pre-cooling the Cold Compress Pack

Daily Storage: After cleaning and drying, store the cold compress pack in the refrigerator compartment (i.e., the "normal temperature compartment"). **Do not store at room temperature outside the refrigerator**, as this may cause the internal coolant to evaporate quickly and draw in air, creating bubbles in the cold compress pack and affecting its performance.

Precautions

Avoid contact with sharp objects to prevent damage and leakage. Do not wash with boiling water or heat in a microwave, as the PVC material is not heat-resistant. Regularly check the seal of the cold compress pack. If damage is found, stop using immediately and replace it.

Disposal Guidelines

Do not dispose of any part of this product as ordinary household waste. This product and its remote control contain lithium-ion polymer batteries. To protect the environment, at the end of the product's service life: Return the fully packaged device to the manufacturer or its authorized representative, or strictly follow the electronic waste disposal regulations in your area for safe recycling.

Common Problems and Solutions

Problem	Solution
The remote control screen does not light up, and the operation is unresponsive, unable to control the mask	- Check if the remote control battery is depleted; use the charging cable to connect to the Type-C port for charging. - If it still does not work after charging, contact the manufacturer or after-sales service center for inspection.
The remote control screen does not light up, but it can control the mask	Reconnect the remote control and the mask following the pairing steps - Press and hold the mask power button for 5 seconds. The indicator will first show a green light, then flash red, indicating the device has entered pairing mode. - Press the power button on the remote control to complete the pairing.
The product cannot be charged	- Replace the charging cable. If the issue is confirmed to be with the charging cable, contact us for a replacement. If the problem persists, contact us to replace the product. Check the charging adapter and voltage. Please ensure that the charging adapter connected to the device is functioning properly and that the voltage specified on the product matches the local mains voltage.
Air bubbles in the cold compress pack / 50% reduction in cold compress pack weight	- Replace with a new cold compress pack

EN Product Parameters

Product Name	Dreame Chrona Pro Photon Anti-Aging Mask
Product Model	AB01P
Battery Parameters	Type: 3.7V / 1000mAh Lithium Polymer Battery (Each) Quantity: 2 Removable Battery Packs (Total Capacity: 2000mAh)
Product Dimensions	Mask: 634*233.8*5.5mm Remote Control: 133.6*42.3*19mm
Net Weight	Mask: Approx. 340.5g Remote Control: Approx. 72g Full-face cold compress pack: approx. 155g Local cold compress pack: approx. 18g
Charging time	Approx. 2 hours
Battery life	Single fully charged battery: approx. 60 minutes Two batteries used together: approx. 120 minutes total

Product Name	Dreame Chrona Pro Photon Anti-Aging Neck Device
Product Model	AB02A
Battery Parameters	Type: 3.7V / 1000mAh built-in lithium polymer battery
Product Dimensions	Neck device: 634*233.8*5.5mm
Net Weight	Neck device: approx. 170g
Charging time	Approx. 2 hours
Battery life	Approx. 60 minutes

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Manufacturer: Dreame Trading (Tianjin) Co., Ltd.