# Oclean

Product User Manual

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
# Oclean Sonic Electric Toothbrush User Manual


⑦

⑧

⑨

⑩

You can read the component description corresponding to each label in detail.

## Thanks for using Oclean Sonic Electric Toothbrush.

**Important Safety Notice**

• When using the product for the first time, due to the high-speed swinging of brush, if it feels stimulating or numb, please do not worry, it is normal and the feeling will weaken and disappear gradually in around 7-14 days.
• If not using for long time (the product battery has no power due to the standby power consumption), when charging for the toothbrush again, the charging indicator doesn't show, please, wait for 6-10 minutes patiently, it can recover normal charging.
• If the main part of toothbrush becomes stony hot when charging, it is normal. If you feel toothbrush to Oclean's After-sales Services Center for testing.
• The product can only use the spect brush head for Oclean Sonic electric toothbrush. Oclean will not be held liable for all the consequences caused by the use of other unofficial or non-licensed brush head products.
• The product is only used for cleaning teeth and must not be used for other purposes.
• The product is only intended as a personal cleaning product. It is not recommended to use it by multiple persons simultaneously.
• The procuct is only used indoors, do not place the toothbrush under the blazing sun for quite a long time. It shall be placed in a dry and well-ventilated place with normal temperature.

**Attention should be paid for those with the following conditions of oral sensitivities**

• Those who have recently undergone periodontal surgery or are undergoing treatments.
• Individuals with other symptoms in the mouth including those with unstable dentures or have persistent heavy bleeding from the gums.
• Other individuals with sensitive and fragile oral issues. This product is not recommended for use with the above group or consult a professional dentist before use.
*This product is not suitable for people with disabilities, mental illness or lack of relevant experience and knowledge [including children (0-12 years old)]. Unless used under the supervision of their guardians.

**Instructions for use**

① ON/OFF button control:
Enable: short press the button
Switch mode: Within 5 seconds after the button is short pressed, you can switch between modes by rapidly short pressing the button when the mode indicator is in the state of "breathing light'.
Turn off:
1. You can turn this product off by short pressing the button again after circularly short pressing the button until this product operates in the "soft mode".
2. You can tum off this product by short pressing the button again when the mode indicator isin the 'normally on' state within 5 seconds after short pressing the
button to start up the toothbrush. '
3. This product is automatically turned off after it is in operation for 2 minutes after you short press the button to start up it. Low battery reminder: When the battery level is lower than 10%, the button circular lamp is in the state of "flicker-ing in red' after you short press the button to turn on the toothbrush.

**Default brushing program:**

Five toothbrushing modes are built in this product: morning mode, night mode, cleaning mode, whitening mode and soft mode.

② Morning mode: compared with “night mode’, the pressure is low. AS a cleaning and care tool, it can effectively remove dental plaque caused during the sleep.
③ Night Mode: deeply remove food residues and dental plaque on the teeth with high pressure.
④ Cleaning mode: meet requirements for routine and ordinary cleaning. Compared with “night mode’, the pressure is relatively tender.
⑤ Whitening mode: effectively clean bite marks in those who have habits of smoking or frequently drink coffee, tea, etc.
⑥ Soft mode: protect the health of gingiva while fulfilling the fundamental function of removing dental plaque in those who have oral allergy.

**Type-C direct connection charging**
Open the bottom cover of the electric toothbrush, connect Type-C interface into the port at the bottom of the toothbrush, connect one end of the USB interface to the power supply, the button circular lamp flickers in white and this indicates that the toothbrush is being charged.

**Charging process**
Charging status: the button circular lamp is on in white and flickers.

**Use for the first time**
The user may charge the toothbrush depending on his or her need for battery life to ensure normal use of this product.

**Application process**
When the battery level is lower than 10%, the button circular lamp is in the state of "flickering in red" after you short press the button to turn on the toothbrush.

**Precautions**
Battery can not be replaced by users themselves, it can only be replaced by the personnel in the customer center autho-rized by Oclean, non-specialized personnel can not dismantle the product by themselves, or it will be regarded as that the user give up the warranty automatically. Please charge the Oclean toothbrush with a 5V/1A or 5V/2A mobile phone charger. Such mobile phone charger includes iPhone Ses tablet charger, etc. If the charging current is lower than 1A, this may cause the charging time to be extended; if you charge the toothbrush with a mobile phone charger with a voltage of higher than 5V, this may cause the charger circuit to be damaged.

**Recommended tooth brushing method**
⑦ Move the toothbrush slowly at 45° between the brush and gingiva line. There is no need to manually brush up and down during this process.
⑧ Hold the toothbrush handle gently, move slowly in a circular motion, and please do not apply excessive pressure.
⑨ Please gargle with clear water after brushing.

**Cleaning and maintenance**
⑩ Rinse with tap water and air dry, and do not pull the bristles of the brush head.

⚠

1. Do NOT disassemble, strike, extrude or throw the device into a fire.
2. Do NOT use the toothbrush if its battery is bulging or has expanded.
3. Do NOT subject the battery to a high-temperature environment.
4. Do NOT use the battery if it has n subjected to moisture.

Before the appliance is discarded, its battery must be removed from it Before the battery is removed, the appliance must be powered off The battery shall be disposed of safely

## • Specifications

| | | | |
|---|---|---|---|
| Product Name | Oclean Sonic Electric Toothbrush | Input | 5B ⎓ 0.8A |
| Model | F5002 | Service Time* | <180 Days (*) |
| IP Rating | IPX7 | Battery Capacity | 2000 mAh(**) |
| Power | 4w | Charging Time | < 5 hours(***) |
| Charger | Type-C Direct Recharge | Implementation Standard | GB39669-2020 |

* Validated by 2 brushing processes per day and 2 minutes of each process. (Morning mode)
** The battery capacity theoretically has a capacity tolerance of +50mAh.
*** It is suggested to adopt the USB charger of 5V = 1A in the phone industry.

Manufacturer:Shenzhen Yunding Information Technology Go.ltd.
Manufacturer's address; Room 412, Building 4, Shenzhen Software Industry Base, 17-19 Haitian Ist Road, Binhai Community, Yuehai Subdistrict, Nanshan District, Shenzhen, Guangdong, China

**FCC Warning**
This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Note: This equipment has been tested and found to comply with the limits for a Class B Satal device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee the interference will not occur in a particular installation. if this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:
— Reorient or relocate the receiving antenna.

— Increase the separation between the equipment and receiver.

— Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
— Consult the dealer or an experienced radio/TV technician for help.

Any changes or modifications not expressly approved by the party responsible for compliance could void your authority to operate the equipment.
When removing the battery, the product must be powered off. After removing the battery, it should be safely disposed of. Please have professional personnel follow the steps below to remove the battery from the product.

**CE Warning**

1. the appliance is not to be used by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction
2. children being supervised not to play with the appliance
3. Signal words WARNING, CAUTION, DANGER in uppercase having a height as specified: Uppercase letter of the text explaining the signal word not smaller than 1,6 mm; Moulded in, engraved, or stamped markings either raised above or have a depth below the surface of at least 0,25 mm, unless
4. This appliance contains batteries that are non-replaceable.

When removing the battery, the product must be powered off. After removing the battery, it should be safely disposed of. Please have professional personnel follow the steps below to remove the battery from the product.

Before abandoning the product, the battery must be removed from the product according to the following steps:

1. First blow the charging board with a Heat gun, and then pry the charging board along the edge of the bottom cover with tweezers.
2. Use a Phillips screwdriver to remove the fixing screws from the bottom cover and remove the bottom cover.
3. Use a nickel to pry the clips on both sides inward and forcefully push the motor shaft until it comes out of the movement.
4. Use scissors to cut off the charging cable at the bottom of the battery, break open the clip, and remove the battery.

– the battery must be removed from the appliance before it is scrapped;
–the appliance must be disconnected from the supply mains when removing the battery;
–the battery is to be disposed of safely.

CAN ICES-003(B)/NMB-003(B)
FDA Registration Number: 3014246581
Made in China

## Warranty Certificate of Oclean Sonic Electric Toothbrush

**User Information**

Name: _____ Contact number: _____ E-mail: _____

ZIP code: _____ Address: _____

_____

**Retailer Information**

Retailer: _____ Sales date: _____

Phone: _____ ZIP code: _____

Address: _____ Invoice no: _____

**Product Information**

Product type: _____

Manufacturing site: _____

**Warranty Card of Oclean Smart Sonic Electric Toothbrush**

| | |
|---|---|
| Delivery date | |
| Fault description/reason/disposal | |
| Equipment employed | |
| Acceptance date | |
| Serviceman | |

Hereby, Shenzhen Yunding Information Technology Co.,Ltd. declares that the radio equipment type F5002 is in compliance with Directive 2014/30/EU, 2014/35/EU, 2011/65/EU. The full text of the EU declaration of conformity is available at the following internet address: http://www.oclean.com/en/service.aspx Please request a maintenance report from the maintenance center and keep it together with the warranty certificate. Shenzhen Yunding Information Technology Co., Ltd., in strict accordance with the Responsibility Regulations on Maintenance, Replacement and Sales Return of Household Appliances, offers a guarantee for consumers who purchase the Oclean Sonic Electric Toothbrush. In the case of any inconsistency between the service articles and the provisions of guarantee, the provisions of the guarantee and the contents of the guarantee certificate shall govern. If the product you have purchased is found to have a fault or defect, please contact our service line at service@oclean.com Website: www.oclean.com